



## 3Shape TRIOS®

### Connect Your 3Shape TRIOS® to DDL's Workflow

Save time and money when you send scans directly from your 3Shape Trios® digital impression scanner.  
**In addition to all-ceramic restorations, DDL accepts digital impressions for PFM and full-cast restorations.**

- In your 3Shape Communicate Account, add 3shape@diversedentallab.com to your list of labs.
- Diverse Dental Lab must authorize your connection to our workflow before you can transfer files.
- Once we confirm your request, visit Laboratory Selection and select Diverse Dental Lab as your laboratory.

If you have any questions, contact DDL at 1.844.463.4837 and ask for digital impression technical support.

### Looking for More Information?

For product information or resources visit [www.diversedentallab.com](http://www.diversedentallab.com) or **Call us on 1.844.463.4837**

Save money and time on every digital case, choose to go model free!

**SAVE 20%**  
ON YOUR  
FIRST DIGITAL CASE

PROMO CODE:

**DDLDIGITAL**

Digital submissions are simple when you choose Diverse Dental Lab.  
**Write promo code: DDLDIGITAL** on your Rx when you upload your first case and save 20% on your first digital submission.