



iTero® Element

Connect Your iTero® to DDL's Workflow

Save time and money when you send scans directly from your iTero® digital impression scanner.
In addition to all-ceramic restorations, DDL accepts digital impressions for PFM and full-cast restorations.

- Call iTero® tech support at 800-577-8767, select option 2 and ask them to add Diverse Dental Lab as a lab partner.
- Select Diverse Dental Lab as your lab partner and submit your digital file

If you have any questions, contact DDL at 1.844.463.4837 and ask for digital impression technical support.

Looking for More Information?

For product information or resources visit www.diversedentallab.com or **Call us on 1.844.463.4837**

Save money and time on every digital case, choose to go model free!

SAVE 20%
ON YOUR
FIRST DIGITAL CASE

PROMO CODE:

DDL DIGITAL

Digital submissions are simple when you choose Diverse Dental Lab.
Write promo code: DDL DIGITAL on your Rx when you upload your first case and save 20% on your first digital submission.