



MidMark™ True Definition™

Connect Your MidMark™ True Definition™ to DDL's Workflow

Save time and money when you send scans directly from your True Definition™ digital impression scanner. **In addition to all-ceramic restorations, DDL accepts digital impressions for PFM and full-cast restorations.**

- Call MidMark™ by dialing 844-856-1231. Select option 3, and then option 1.
- Request that they add Diverse Dental Lab to your account.
- Log out then log back in to your MidMark™ True Definition™ account.
- Select Diverse Dental Lab from the laboratory options.

If you have any questions, contact DDL at 1.844.463.4837 and ask for digital impression technical support.

Looking for More Information?

For product information or resources
visit www.diversedentallab.com or
Call us on 1.844.463.4837

Save money and time on every digital case,
choose to go model free!

SAVE 20%
ON YOUR
FIRST DIGITAL CASE

PROMO CODE:

DDLDIGITAL

Digital submissions are simple when
you choose Diverse Dental Lab.
Write promo code: DDLDIGITAL on
your Rx when you upload your first case
and save 20% on your first digital
submission.